HIGH THROUGHPUT SCREENING (HTS)

By Dr. Debasish Pradhan

Sr. Faculty Pharmacology

HIGH THROUGHPUT SCREENING (HTS) is identification of one of more positive candidates extracted from a pool of possible candidate based on specific criteria.

LEAD compound is a representative of a compound series with sufficient potential (as measured by potency, selectivity, pharmacokinetic physicochemical properties, toxicity and novelty) to progress to a fundament development programme.

- *It is a drug-discovery process widely used in the pharmaceutics industry.
- *It allows automation to quickly assay the biological or biochemica activity of a large no of compounds.
- *HTS is process by which large nos. of compounds are rapidly tested for their ability to modify the properties of a selected biological target.

It is a useful for discovering ligands for receptors, enzymes, ionchannels or other pharmacological targets, or pharmacologically profiling a cellular or biochemical pathway of interest.

Screening Collection

Screening libraries:

- (a) Diversity Libraries
- (b) Focused Libraries:

Ligand-Based Approaches Structure-Based Approaches

Compound acquisition Purchasing commercially available

LIBRARIES

- *siRNA libraries: providing the best available RNA interference technology with maximum flexibility The mouse whole genome contain about 16872 gene targets and the Human whole-genome contains 18120 gene targets.
- microRNA libraries: Human libraries of synthetic microRNA mimics and microRNA inhibitors.
- *Small Compound libraries: A library of 640 FDA approved drugs (Screen-Well FDA Approved Drug Library, Enzo Life

Steps in HTS

***** 1 st stage screening

- Test optical clarity, abrasion resistance, and adhesion
- Eliminates ~ 90% of samples

2 nd stage screening

- Test weather ability, integrity, gloss, and surface smoothness
- ~10% of the samples

DETECTION METHODS IN HTS

- Spectroscopy
- Mass Spectrometry
- Chromatography
- Calorimetry
- X-ray diffraction
- Microscopy
- Radioactive methods

SPECTROSCOPY IN HTS:

- Fluorescence Spectroscopy
- Total internal reflection fluorescence (TIRF)
- Nuclear magnetic resonance (NMR)
- Absorption and luminescence
- Fourier transformed infrared(FTIR)
- Light scattering

CHROMATOGRAPHY IN HTS:

- Gas chromatography (GC)
- Thin layer chromatography
- Liquid chromatography (HPLC)

CALORIMETRY IN HTS:

- Isothermal Titration Calorimetry (ITC)
- Differential Scanning Calorimetry (DSC)

MICROSCOPY IN HTS:

- Scanning Tunnelling Microscopy
- Atomic Force Microscopy
- Confocal Microscopy

USES:

To screen Micro arrays such as:

- DNA chips
- · RNA chips

METHODOLOGY

The heart of the HTS system is a plate, or tray, which consists of tiny wells where assay reagents and samples are deposited, and their reactions monitored

-The configuration of the plate has changed from 96 wells (in a matrix of 8 rows by 12 columns) to 384, and now to a high - density 1536 - well format, which enables large - scale screening.

-Assay reagents may be coated onto the plates or deposited in liquid form together with test samples into the wells.


Both samples and assay reagents may be incubated and

Cell - based assays have become an important test compared with other *in vitro assays*, as they can provide information about bioavailability, cytotoxicity and effects on biochemical pathway

The enzyme - based and cell - based assay systems consist of receptors or mimetics of receptors (components that mimic active parts of receptors)

Receptors- Normally the assays are linked to an indicator that shows the ligands – receptor interaction as some form of signal

The advantage of cell - based assays over biochemical assays is that cell - based assays enable the analysis of sample compound



CELL BASED ASSAYS

- Cell-based assays refer to any of a number of different experiments based on the use of live cells
- This is a general definition and can include a variety of assays that measure cell proliferation, toxicity, motility, production of a measurable product and morphology
- Cell-based assays offer a more accurate representation of the

FOUR KEY ELEMENTS OF CELL BASED ASSAY:

- A cellular component e.g. a cell line or a primary cell population
- A target (substrate) molecule that records the cellular response
- An instrument to conduct and monitor the assay An informatics component to manage and analyse data from the assay
- Cell-based reporter assays are used where human receptors are transfected into null cell lines either alone, (luciferinluciferase) or light transmission (melanophore), that can be measured independently of radioactivity within minutes.

ADVANTAGES

- Assays do not require purification of the target protein
- Can immediately select against compounds / potential drugs that are generally cytotoxic, or that cannot permeate cellular membranes to reach intracellular sites
- * Hit/lead compounds identified by cell based assays have passed important validation steps, saving time and costs in drug development
- Cell-based assays visualize all possible drug-target interactions e.g. activators, target interactions
- * High sensitivity of assays & Automation

Newer method in HTS

- * High-Throughput, Fluorescence-Based Screening
- (a) For screening, expressed proteins are labelled either as fusions with green fluorescent protein (GFP) or through translational incorporation of a fluorescent amino acid derivative, BODIPY-FL-Lysine.
- (a) Using fluorescence detection, the entire procedure can be carried out in approximately 8 h.

Applications

- In drug discovery
- Systematic Study of Mitochandrial Toxicity of

CELLULAR COMPONENTS

- Different cell lines are being used in cell based assays
- Some examples are:
 - **HUMAN CELL LINES**
 - DU145, PC3, Lncap (Prostate cancer)
 - MCF-7, MDA-MB-438, T47D (Breast cancer)

THD 1 (A and a Manifeld I and and a

REFERENCES

Leach AR, Hann MM. The in silico world of virtual libraries. Drug Discov Today 2000; 5:326–336.

J. AM. CHEM. SOC. 2010, 132, 13182-13184

THANK YOU